


MECCA CATHE

1. DISPOSAL MEDICAL DEVICES


- Central Venous Catheter made of specially formulated and biocompatible Polyurethane material provides strength during insertion and also softens at body temperature to conform to the body tissues and reduces the risk and vascular trauma.
- Specially designed soft & beveled tip for smooth & easy insertion of catheter
- Soft Flexible J-Tip Guide wire prevents the vessel perforation and also provides good torque to ensure film insertion
- Sufficiently radio-Opaque material of catheter with clear, definite marking facilitates correct placement of cathetter tip
- Sterile ready for use
- Individually tray packed

FOUR LUMEN




Specially Designed Soft & Conical Tip ensures Atraumatic, Safe & Easy Insertion

TRIPLE LUMEN



Biocompatible PU material is stiff & non kinking at room temperature, but softens at body temperature

DOUBLE LUMEN




Radiopaque catheter & length markings ensures accurate placement.

SINGLE LUMEN



Full range of Adult, Pediatric & Infant CVC Kits are available.

2. COMPONENTS OF STANDARD KIT



Kink resistant Nickel-Titanium Guidewire with soft & flexible J-tip offers better torque which helps in easy insertion & prevents vessel perforation.
Wires are Double Distal.

2 piece design of guidewire advancer.




CVC kits are available with options of std. Needle, Needle with Back Flow valve & Introducer Syringe.



COMPLETE SET OF CVC KIT CONSIST OF:


- Indwelling catheter
- Y-Introducer needle with check valve
- J-Tip guide wire
- Vessel dilator
- Luer lock syring
- Scalpel
- Catheter holder
- Catheter holder clamp

- Injection cap
- Extension Line Clamp
- Guiding Syringe/LOR Syringe (Optional)



ANGI

1 MANIFOLD



- Choice of 2, 3 or 4 Core
- Options of colour handle
- Clear material for easy detection of bubbles
- Large bore for maximum flow 500 PSI
- Variants available on request.

3. CONTROL SYRINGE



- Choice of 10cc, 12cc or 20cc capacity
- Low friction hub
- Bold marking for correct measurement of dosage
- Finger grips for one hand operation
- Rotating male luer lock

4. CONTRAST LINE



Customised sets available upon order

5. ORDERING INFORMATION



- Choice of length from 10", 20" or 48", 60"
- Braided PU material, 1200 PSI rating
- Rotating male luer lock
- Also available in PVC material upon request

6. GUIDE WIRE



- Variety of guide wires available in following variants:
- O.D. : 0.018", 0.021", 0.025", 0.032", 0.035", 0.038"
- Length : 45cm, 60cm, 70cm, 150cm, 260cm
- Tip : J tip, Straight
- Material : Stainless steel wires with PTFE coated and Non coated, Nickle titanium

7. INTRODUCER NEEDLE



Highly precise, atraumatic bevel tip needle

8. HEMOSTASIS: PUSH-PULL/ ROTATING/ CLICK TYPE



Large bore Compatible with catheters from 4Fr upto 9Fr and guide wires from 0.014" to 0.038" dia Rated pressure 480 PSI
Spcl. feature of Click Y Septum can be used in OPEN, SEMI OPEN and FULL CLOSED position

CUSTOMISED EXTENSION LENGTHS AVAILABLE UPON REQUEST

a) TORQUE DEVICE



Perfect grip for easy manipulation of guidewire

b) INTERSECTION TOOL/ BLUNT NEEDLE



CATHLAB KITS



ORDERING INFORMATION/CONTENTS

- MEC-01 (3RN) Manifold 3 Port R/On, IV Set, Contrast Line, 10ml Control Syringe, 40" PM Line, 'Rotating' Hemostasis Y, Torque Device, Blunt Needle, 20" M/F extension & Injector Connector.
- MEC-02(2RN) Manifold 2 Port R/On, IV Set, 10ml LL Syringe, 40" & 20" PM Line, 10mlsyringes
- MEC-03 Hemostasis Y 'Click-Press' type with 10" extension, Torque Device & Blunt Needle.
- MEC-04 Hemostasis Y, 'Rotating' type, Torque Device, & Blunt Needle
- MEC-05 Accessories Kit Includes Hemostasis Y Click Type, Torque Device, Needle.
- MEC-06 Hemostasis Y, 'Push-Pull' type with 10" extension, Torque Device, & Blunt Needle
- MEC-07 Inflation Device 30atm, Hemostasis Y 'Rotating' type, Torque Device, & Blunt Needle
- MEC-08 Accessories Kit includes Push/Pull Y with Extension 25 cms and 4 way Stopcock.

(Other Customised kits also available upon request)

All the components are available individually packed or in assembled kit.