



MECCA HEALTHCARE

Haemodialysis Catheter Kit

Innovative Technology for Haemodialysis Catheter

Unique design with Soft Tip, Side Holes, Rotatable Suture Wings and Silicone Extension Tube

CARE COMPASSION TRUST


www.mhplindia.com

+91 999 869 9090

MECCA HEALTHCARE HAEMODIALYSIS CATHETER KIT

Features and benefits

- The Catheter is made from medical grade PU, which has better flexibility and anti-kinking performance
- Flexibility & Floppy
- Thermosensitive, softens as it warms to body temperature
- Smooth inside and outside, avoids clot formation
- Radiopaque, facilitates quick visualization under X-ray
- No scale, all of the tube shall be inserted completely



Standard Kit includes:

- Haemodialysis Catheter
- Vessel Dilator
- Introducer Needle
- Syringe
- Guide-wire
- Adhesive Wound Dressings
- Heparin Caps
- Scalpel
- Needle with Suture


Optional Compound Kit includes:

- Standard Kit accessories
- 5ml Syringe
- Surgical Gloves
- Surgical Pledget
- Surgery Sheet

lumen					
	Single				
	Double I				
	Double II				
	Triple	Silicone extension tube	Chronic	Compound kits	Nitinol guide-wire

Soft Tip

Tapered tip facilitates ease of insertion without use of peel-away sheath. Reduces vessel trauma during insertion




Side Holes

2.5 cm far from end hole reduces recircling and wall suction. Minimizes risk of clot formation.




Rotatable Suture Wing

Facilitates skin inspection. Avoid blood vessel wall-suction. Minimizes risk of exit infections. No tube-kinking occurred



Silicone Extension Tube

Increases patient comfort, will not crimp over time, able to visualize fluids.



Our Complete Product list

Anesthesia

P.M.KIT
C.V.P.MANOMETER
OXYGEN MASK
TWIN OXYGEN SET
OXYGEN RECOVERY KIT
VENTILATOR CIRCUIT
CVC KIT
H.M.E.FILTER
CATHETER MOUNT
P.M.LINE
EXTENSION LINE
3WAY STOPCOCK WITH EXT. LINE
3-WAY STOP COCK
BAIN CIRCUIT
E.T.TUBE
SPINAL NEEDLE
GUIDEL AIRWAY

Cardio thoracic and Vascular surgery

HEART LUNG BYPASS TUBING SET
ONLINE CARDIOPLEGIA DELIVERY SYSTEM
ARTERIAL LINE FILTER ADULT
A.V.TUBING
CARDIAC CONNECTORS
CHEST DRAINAGE BOTTLE
C.D.CATHETERS
SUCKER
CARDIOPLEGIA ADAPTER

General Surgery

FOLEY CATHETER
YANKAUER HANDLE
CLOSED WOUND DRAINAGE SET
ARTERIAL DRAINAGE KIT

General Medical Disposables

I.V.SET
B.T.SET
M.V.SET
SUCTION CATHETER
RYLE'S TUBE
INFANT FEEDING TUBE
MUCUS EXTRACTOR
URINE COLLECTION BAG
I.V.CANNULA
LUNG EXERCISER

Cardiology

ANGIO KIT
PTCA KIT
MANIFOLDS
INTRODUCER NEEDLE
INTRODUCER SET
CONTROL SYRINGE
GUIDE WIRE

Gynaecology

M.T.P.CANNULA
UMBILICAL CATHETER
UMBILICAL CORD CLAMP

Urology

TUR SET
MALE CATHETER
NELATON CATHETER
URETHRAL CATHETER

Dialysis

P.D.SET
P.D.CATHETER
HAEMODIALYSIS CATHETER



Mecca Care

A Technical Collaboration with ACME UK Inc. Limited London, Great Britain

Marketed by Mecca Healthcare Pvt. Limited

102, GIDC, Kalol (N.G.), Gandhinagar, Gujarat, India.

For Demonstration purpose contact: 079 4800 5989
www.mhplindia.com