




Meca Care



Nephrology Products

MECCA HEALTHCARE P. LIMITED



Guidewires

Guidewires


We the MECCA HEALTHCARE P LIMITED, as manufacturer

We have worked hard towards achieving excellence in the process of quality "GUIDE WIRES in which we have undergone 100,000 hours of skilled manpower. In the year 2012, we got ourselves collaborated with market leader from Ireland. Sophisticated Technology, Quality Machineries and Expertise Research have been inculcated in our processing unit in Kalol, India. This venture help us to accomplish our goal to make remarkable products" Guide Wire, the metal and raw materials are being imported from renowned companies situated in countries like Ireland, Netherlands and United State of America.

MECCA HEALTHCARE P LIMITED are Indian pioneers to introduce Coated Guidewires, Uncoated Guidewires, Nitinol Guidewires, Hydrophilic Guidewires and Lunderquist Guidewires. Our Guidewires are available with J tip or straight tip or double distal ends. We have SS Core Stiff Guidewires, Extra Stiff Guidewires and Super Stiff Guidewires to meet every Clinical requirement of Professionals. We do have Kink Resistive guidewires as well.

Our Guidewires range is available from diameter 0.018" to 0.038", which supports length up-to 260 cm.

Product Feature



Dialysis Guidewires : Standard

REF	Specification	Coating
MI/GW/3270JUN	0.032 Inch ,70cm, J Tip	Uncoated *
MI/GW/3270STUN	0.032 Inch ,70cm, ST Tip	Uncoated *
MI/GW/3570JUN	0.035 Inch ,70cm, J Tip	Uncoated *
MI/GW/3570STUN	0.035 Inch ,70cm, ST Tip	Uncoated **
MI/GW/3870JUN	0.038 Inch ,70cm, J Tip	Uncoated **
MI/GW/3870STUN	0.038 Inch ,70cm, ST Tip	Uncoated *
MI/GW/3280JUN	0.032 Inch ,80cm, J Tip	Uncoated
MI/GW/3280STUN	0.032 Inch ,80cm, ST Tip	Uncoated
MI/GW/3580JUN	0.035 Inch ,80cm, J Tip	Uncoated
MI/GW/3580STUN	0.035 Inch ,80cm, ST Tip	Uncoated
MI/GW/3880JUN	0.038 Inch ,80cm, J Tip	Uncoated
MI/GW/3880STUN	0.038 Inch ,80cm, ST Tip	Uncoated
MI/GW/3270J	0.032 Inch ,70cm, J Tip	PTFE Coating
MI/GW/3270ST	0.032 Inch ,70cm, ST Tip	PTFE Coating
MI/GW/3570J	0.035 Inch ,70cm, J Tip	PTFE Coating
MI/GW/3570ST	0.035 Inch ,70cm, ST Tip	PTFE Coating
MI/GW/3870J	0.038 Inch ,70cm, J Tip	PTFE Coating
MI/GW/3870ST	0.038 Inch ,70cm, ST Tip	PTFE Coating
MI/GW/3280J	0.032 Inch ,80cm, J Tip	PTFE Coating
MI/GW/3280ST	0.032 Inch ,80cm, ST Tip	PTFE Coating
MI/GW/3580J	0.035 Inch ,80cm, J Tip	PTFE Coating
MI/GW/3580ST	0.035 Inch ,80cm, ST Tip	PTFE Coating
MI/GW/3880J	0.038 Inch ,80cm, J Tip	PTFE Coating
MI/GW/3880ST	0.038 Inch ,80cm, ST Tip	PTFE Coating

*Standard Configuration : 3mm St or J Tip, Standard Stiffness, Fixed Core, PTFE Coating over SS.

Dialysis Guidewires : Kink Resistive

REF	Specification	Coating
MI/GWN/3270JUN	0.032 Inch ,70cm, J Tip	Uncoated
MI/GWN/3270STUN	0.032 Inch ,70cm, ST Tip	Uncoated
MI/GWN/3570JUN	0.035 Inch ,70cm, J Tip	Uncoated
MI/GWN/3570STUN	0.035 Inch ,70cm, ST Tip	Uncoated
MI/GWN/3870JUN	0.038 Inch ,70cm, J Tip	Uncoated
MI/GWN/3870STUN	0.038 Inch ,70cm, ST Tip	Uncoated
MI/GWN/3280JUN	0.032 Inch ,80cm, J Tip	Uncoated
MI/GWN/3280STUN	0.032 Inch ,80cm, ST Tip	Uncoated
MI/GWN/3580JUN	0.035 Inch ,80cm, J Tip	Uncoated
MI/GWN/3580STUN	0.035 Inch ,80cm, ST Tip	Uncoated
MI/GWN/3880JUN	0.038 Inch ,80cm, J Tip	Uncoated
MI/GWN/3880STUN	0.038 Inch ,80cm, ST Tip	Uncoated
MI/GWN/3270J	0.032 Inch ,70cm, J Tip	PTFE Coating
MI/GWN/3270ST	0.032 Inch ,70cm, ST Tip	PTFE Coating
MI/GWN/3570J	0.035 Inch ,70cm, J Tip	PTFE Coating
MI/GWN/3570ST	0.035 Inch ,70cm, ST Tip	PTFE Coating
MI/GWN/3870J	0.038 Inch ,70cm, J Tip	PTFE Coating
MI/GWN/3870ST	0.038 Inch ,70cm, ST Tip	PTFE Coating
MI/GWN/3280J	0.032 Inch ,80cm, J Tip	PTFE Coating
MI/GWN/3280ST	0.032 Inch ,80cm, ST Tip	PTFE Coating
MI/GWN/3580J	0.035 Inch ,80cm, J Tip	PTFE Coating
MI/GWN/3580ST	0.035 Inch ,80cm, ST Tip	PTFE Coating
MI/GWN/3880J	0.038 Inch ,80cm, J Tip	PTFE Coating
MI/GWN/3880ST	0.038 Inch ,80cm, ST Tip	PTFE Coating

*Standard Configuration : 3mm St or J Tip, Standard Stiffness, Fixed Core, PTFE Coating over SS.

Angiography Guidewires : Standard

REF	Specification	Coating
MI/GW/18150J	0.018 Inch ,150cm, J Tip	PTFE Coating *
MI/GW/18150ST	0.018 Inch ,150cm, St Tip	PTFE Coating *
MI/GW/25150J	0.025 Inch ,150cm, J Tip	PTFE Coating *
MI/GW/25150ST	0.025 Inch ,150cm, St Tip	PTFE Coating *
MI/GW/28150J	0.028 Inch ,150cm, J Tip	PTFE Coating *
MI/GW/28150ST	0.028 Inch ,150cm, St Tip	PTFE Coating *
MI/GW/32150J	0.032 Inch ,150cm, J Tip	PTFE Coating **
MI/GW/32150ST	0.032 Inch ,150cm, St Tip	PTFE Coating **
MI/GW/35150J	0.035 Inch ,150cm, J Tip	PTFE Coating **
MI/GW/35150ST	0.035 Inch ,150cm, St Tip	PTFE Coating **
MI/GW/38150J	0.038 Inch ,150cm, J Tip	PTFE Coating **
MI/GW/38150ST	0.038 Inch ,150cm, St Tip	PTFE Coating **
MI/GW/18260J	0.018 Inch ,260cm, J Tip	PTFE Coating *
MI/GW/18260ST	0.018 Inch ,260cm, St Tip	PTFE Coating *
MI/GW/25260J	0.025 Inch ,260cm, J Tip	PTFE Coating *
MI/GW/25260ST	0.025 Inch ,260cm, St Tip	PTFE Coating *
MI/GW/28260J	0.028 Inch ,260cm, J Tip	PTFE Coating *
MI/GW/28260ST	0.028 Inch ,260cm, St Tip	PTFE Coating *
MI/GW/32260J	0.032 Inch ,260cm, J Tip	PTFE Coating *
MI/GW/32260ST	0.032 Inch ,260cm, St Tip	PTFE Coating *
MI/GW/35260J	0.035 Inch ,260cm, J Tip	PTFE Coating *
MI/GW/35260ST	0.035 Inch ,260cm, St Tip	PTFE Coating *
MI/GW/38260J	0.038 Inch ,260cm, J Tip	PTFE Coating *
MI/GW/38260ST	0.038 Inch ,260cm, St Tip	PTFE Coating *

*Standard Configuration : 3mm St or J Tip, Standard Stiffness, Fixed Core, PTFE Coating over SS.

No Star - No Sampling, Check Availability before confirming order of Non Running products.

* No Sampling Available, Check Availability of Running Products.

** Sampling Available, No need to check availability of highly moving products, Easy Available.

Hydrophilic Guide Wire

REF	Specification	Coating
MI/GWH/25150J	0.025 Inch ,150cm, J Tip	HYDRO Coating *
MI/GWH/25150ST	0.025 Inch ,150cm, St Tip	HYDRO Coating *
MI/GWH/28150J	0.028 Inch ,150cm, J Tip	HYDRO Coating *
MI/GWH/28150ST	0.028 Inch ,150cm, St Tip	HYDRO Coating *
MI/GWH/32150J	0.032 Inch ,150cm, J Tip	HYDRO Coating *
MI/GWH/32150ST	0.032 Inch ,150cm, St Tip	HYDRO Coating *
MI/GWH/35150J	0.035 Inch ,150cm, J Tip	HYDRO Coating *
MI/GWH/35150ST	0.035 Inch ,150cm, St Tip	HYDRO Coating *
MI/GWH/38150J	0.038 Inch ,150cm, J Tip	HYDRO Coating *
MI/GWH/38150ST	0.038 Inch ,150cm, St Tip	HYDRO Coating *
MI/GWH/18260J	0.018 Inch ,260cm, J Tip	HYDRO Coating *
MI/GWH/18260ST	0.018 Inch ,260cm, St Tip	HYDRO Coating *
MI/GWH/25260J	0.025 Inch ,260cm, J Tip	HYDRO Coating *
MI/GWH/25260ST	0.025 Inch ,260cm, St Tip	HYDRO Coating *
MI/GWH/28260J	0.028 Inch ,260cm, J Tip	HYDRO Coating *
MI/GWH/28260ST	0.028 Inch ,260cm, St Tip	HYDRO Coating *
MI/GWH/32260J	0.032 Inch ,260cm, J Tip	HYDRO Coating *
MI/GWH/32260ST	0.032 Inch ,260cm, St Tip	HYDRO Coating *
MI/GWH/35260J	0.035 Inch ,260cm, J Tip	HYDRO Coating *
MI/GW/35260ST	0.035 Inch ,260cm, St Tip	HYDRO Coating *
MI/GW/38260J	0.038 Inch ,260cm, J Tip	HYDRO Coating *
MI/GW/38260ST	0.038 Inch ,260cm, St Tip	HYDRO Coating *

Short Term Dialysis Catheter Kit



Product Feature



Soft Tip

MHPL Short term dialysis catheter has been able to Reduces trauma to vessel wall upon introduction and in-situ. The radiopaque catheter and tip facilitates vessel location when visualized by X-Ray.



Curved Extensions

The extension are curved for jugular vein approaches to make the catheter more comfortable for the patient.



Biocompatible Polyurethane Material

With new generation Biocompatible polyurethane material which is tested for toxicity and hemo-compatability and this become softer in body.



Rotating Suture Wing

Rotating suture wing provides secure external anchoring catheter/vein interface .



Antimicrobial Heparin Cap

Specially made adaptor helps to avoid infection.



Multi-Lumen Available

The short term dialysis catheter is available single, double, triple lumen for infant and adult.

Specifications

Short Term Dialysis Catheter

REF	Lumen	French	Length	Inner Diameter	Type
MI/HDC/SL/0713	Single	7.0 Fr	13 cm	14G	Straight
MI/HDC/SL/0715	Single	7.0 Fr	15 cm	14G	Straight
MI/HDC/SL/0716	Single	7.0 Fr	16 cm	14G	Straight
MI/HDC/SL/0720	Single	7.0 Fr	20 cm	14G	Straight
MI/HDC/SL/8513	Single	8.5 Fr	13 cm	14G	Straight
MI/HDC/SL/8515	Single	8.5 Fr	15 cm	14G	Straight
MI/HDC/SL/8516	Single	8.5 Fr	16 cm	14G	Straight
MI/HDC/SL/8520	Single	8.5 Fr	20 cm	14G	Straight
MI/HDC/DL/8511	Double	8.5 Fr	11 cm	14G X 14G	Straight
MI/HDC/DL/8511J	Double	8.5 Fr	11 cm	14G X 14G	Curved
MI/HDC/DL/11513	Double	8.5 Fr	13 cm	12G X 12G	Straight
MI/HDC/DL/11513J	Double	11.5 Fr	13 cm	12G X 12G	Curved
MI/HDC/DL/11515	Double	11.5 Fr	15 cm	12G X 12G	Straight
MI/HDC/DL/11515J	Double	11.5 Fr	15 cm	12G X 12G	Curved
MI/HDC/DL/11516	Double	11.5 Fr	16 cm	12G X 12G	Straight
MI/HDC/DL/11516J	Double	11.5 Fr	16 cm	12G X 12G	Curved
MI/HDC/DL/11520	Double	11.5 Fr	20 cm	12G X 12G	Straight
MI/HDC/DL/11520J	Double	11.5 Fr	20 cm	12G X 12G	Curved
MI/HDC/DL/1213	Double	11.5 Fr	13 cm	12G X 12G	Straight
MI/HDC/DL/1213J	Double	12.0 Fr	13.5 cm**	12G X 12G	Curved
MI/HDC/DL/1215	Double	12.0 Fr	15 cm**	12G X 12G	Straight
MI/HDC/DL/1215J	Double	12.0 Fr	15 cm**	12G X 12G	Curved
MI/HDC/DL/1216	Double	12.0 Fr	16 cm**	12G X 12G	Straight
MI/HDC/DL/1216J	Double	12.0 Fr	16 cm**	12G X 12G	Curved
MI/HDC/DL/1220	Double	12.0 Fr	20 cm**	12G X 12G	Straight
MI/HDC/DL/1220J	Double	12.0 Fr	20 cm*	12G X 12G	Curved
MI/HDC/TL/1213	Triple	12.0 Fr	13.5 cm	16G X 12G X 12G	Straight
MI/HDC/TL/1213J	Triple	12.0 Fr	13.5 cm**	16G X 12G X 12G	Curved
MI/HDC/TL/1215	Triple	12.0 Fr	15 cm**	16G X 12G X 12G	Straight
MI/HDC/TL/1215J	Triple	12.0 Fr	15 cm**	16G X 12G X 12G	Curved
MI/HDC/TL/1216	Triple	12.0 Fr	16 cm	16G X 12G X 12G	Straight
MI/HDC/TL/1216J	Triple	12.0 Fr	16 cm	16G X 12G X 12G	Curved
MI/HDC/TL/1220	Triple	12.0 Fr	20 cm	16G X 12G X 12G	Straight
MI/HDC/TL/1220J	Triple	12.0 Fr	20 cm	16G X 12G X 12G	Curved

Short Term Dialysis Catheter Kit

REF	Lumen	French	Length	Inner Diameter	Type
MI/HDCK/SL/0713	Single	7.0 Fr	13 cm	14G	Straight
MI/HDCK/SL/0715	Single	7.0 Fr	15 cm	14G	Straight
MI/HDCK/SL/0716	Single	7.0 Fr	16 cm	14G	Straight
MI/HDCK/SL/0720	Single	7.0 Fr	20 cm	14G	Straight
MI/HDCK/DL/8511	Double	8.5 Fr	11 cm	14G X 14G	Straight
MI/HDCK/DL/8511J	Double	8.5 Fr	11 cm	14G X 14G	Curved
MI/HDCK/DL/11513	Double	11.5 Fr	13 cm	12G X 12G	Straight
MI/HDCK/DL/11513J	Double	11.5 Fr	13 cm	12G X 12G	Curved
MI/HDCK/DL/11515	Double	11.5 Fr	15 cm	12G X 12G	Straight
MI/HDCK/DL/11515J	Double	11.5 Fr	15 cm	12G X 12G	Curved
MI/HDCK/DL/11516	Double	11.5 Fr	16 cm	12G X 12G	Straight
MI/HDCK/DL/11516J	Double	11.5 Fr	16 cm	12G X 12G	Curved
MI/HDCK/DL/11520	Double	11.5 Fr	20 cm	12G X 12G	Straight
MI/HDCK/DL/11520J	Double	11.5 Fr	20 cm	12G X 12G	Curved
MI/HDCK/DL/1213	Double	12.0 Fr	13 cm	12G X 12G	Straight
MI/HDCK/DL/1213J	Double	12.0 Fr	13 cm	12G X 12G	Curved
MI/HDCK/DL/1215	Double	12.0 Fr	15 cm	12G X 12G	Straight
MI/HDCK/DL/1215J	Double	12.0 Fr	15 cm	12G X 12G	Curved
MI/HDCK/DL/1216	Double	12.0 Fr	16 cm	12G X 12G	Straight
MI/HDCK/DL/1216J	Double	12.0 Fr	16 cm	12G X 12G	Curved
MI/HDCK/DL/1220	Double	12.0 Fr	20 cm	12G X 12G	Straight
MI/HDCK/DL/1220J	Double	12.0 Fr	20 cm	12G X 12G	Curved
MI/HDCK/TL/1213	Triple	12.0 Fr	13 cm	16G X 12G X 12G	Straight
MI/HDCK/TL/1213J	Triple	12.0 Fr	13 cm	16G X 12G X 12G	Curved
MI/HDCK/TL/1215	Triple	12.0 Fr	15 cm	16G X 12G X 12G	Straight
MI/HDCK/TL/1215J	Triple	12.0 Fr	15 cm	16G X 12G X 12G	Curved
MI/HDCK/TL/1216	Triple	12.0 Fr	16 cm	16G X 12G X 12G	Straight
MI/HDCK/TL/1216J	Triple	12.0 Fr	16 cm	16G X 12G X 12G	Curved
MI/HDCK/TL/1220	Triple	12.0 Fr	20 cm	16G X 12G X 12G	Straight
MI/HDCK/TL/1220J	Triple	12.0 Fr	20 cm	16G X 12G X 12G	Curved

No Star - No Sampling, Check Availability before confirming order of Non Running products.

* No Sampling Available, Check Availability of Running Products.

** Sampling Available, No need to check availability of highly moving products, Easy Available.



LongTerm Dialysis Catheter

We develop, manufacture, market and support cutting edge vascular access devices and accessories to meet the clinical needs of the medical industry, particularly in the field of interventional cardiology and dialysis. Our company's engineering and applications expertise provides superior products whose progressive designs accommodate advances in medical device and whose quality anticipates the requirements of our professional clients and the patients they serve.

Currently one of the largest manufacturer & importer in the field of Cardiology, Nephrology & Surgical Disposables in INDIA.

LongTerm Dialysis Catheter

MECCA HEALTHCARE P LIMITEDs dialysis Catheter is a reliable long-term dialysis catheter that consistently delivers high flow rates with low venous and arterial pressures.

Standard Kit Contains :

- Longterm Dialysis Catheter
- Vessel Dilator
- Valved Peelable Introducer
- Introducer Needle
- Guide-wire
- Tunneler (Tocar)
- Adhesive Wound Dressings
- Herparin Caps
- Gauze Sponges
- Scalpel
- Hollow Syringe
- Needle with Suture



Double D Design

Long term dialysis catheter is made of double D design tubing, which provide excellent flow with low pressure.


Various of Lengths & tip-to-cuff distance

Our long range of catheters help in securing the desired catheter for diverse patient sizes.



Carbothane Material

With second generation Carbothane polyurethane material which facilitates percutaneous placement also this resists collapse during treatment. Which also gives our catheter excellent flexibility.




Infinity tip transition

Our Peel away catheter comes with unique infinity tip transition, which provide excellent insertion and helps in minimizing trauma .



Cuff

MHPL Long term dialysis catheter has 9mm of cuff length.



New Tunneler

With atraumatic tip and back end our tunneler helps in minimizing trauma.

Specifications

Double Lumen Long term Dialysis Catheter Kit

REF	French	Catheter From Cuff	Length From Hub	A. Volume	V. Volume
MI/PCK/14519ST	14.5 Fr	19 cm	24 cm	1.4 cc	1.8 cc
MI/PCK/14523ST	14.5 Fr	23 cm	28 cm	1.8 cc	2.2 cc

No Star - No Sampling, Check Availability before confirming order of Non Running products.

* No Sampling Available, Check Availability of Running Products.

** Sampling Available, No need to check availability of highly moving products, Easy Available.

We also offer Short Term Dialysis Catheters, Fistulla Needles, Transducer Protector, Dialysis Wire and Introducer Needles.



Single Lumen Femoral Catheter

We develop, manufacture, market and support cutting edge vascular access devices and accessories to meet the clinical needs of the medical industry, particularly in the field of interventional cardiology and dialysis. Our company's engineering and applications expertise provides superior products whose progressive designs accommodate advances in medical device and whose quality anticipates the requirements of our professional clients and the patients they serve.

Currently one of the largest manufacturer & importer in the field of Cardiology, Nephrology & Surgical Disposables in INDIA.


Single Lumen™ Femoral Catheter

MECCA HEALTHCARE P LIMITEDs offers Catheter & Kit. It covers almost every size in pediatric and adult length. It comes in single lumen, double lumen & triple lumen with kink resistive catheter.

Standard Kit Contains :

- Guidewire
- Single Lumen Catheter
- Introducer Needle
- Injection Cap
- Dialator

Product Feature



Double D Design

Long term dialysis catheter is made of double D design tubing, which provide excellent flow with low pressure.




Various of Lengths & tip-to-cuff distance

Our long range of catheters help in securing the desired catheter for diverse patient sizes.



Carbothane Material

With second generation Carbothane polyurethane material which facilitates percutaneous placement also this resists collapse during treatment. Which also gives our catheter excellent flexibility.



Infinity tip transition

Our Peel away catheter comes with unique infinity tip transition, which provide excellent insertion and helps in minimizing trauma .



Cuff

Long term dialysis catheter has 9mm of cuff length.



New Tunneler

With atraumatic tip and back end our tunneler helps in minimizing trauma.

Specifications

Single Lumen Femoral Catheter

REF	Catheter Length Specification
MI/UDC/SL/01	7 Fr, ID:14G Length:13 cm
MI/UDC/SL/02	7 Fr, ID:14G Length:15 cm
MI/UDC/SL/03	7 Fr, ID:14G Length:16 cm
MI/UDC/SL/04	7 Fr, ID:14G Length:20 cm
MI/UDC/SL/05	8 Fr, ID:12G Length:13 cm
MI/UDC/SL/06	8 Fr, ID:12G Length:15 cm
MI/UDC/SL/07	8 Fr, ID:12G Length:16 cm
MI/UDC/SL/08	8 Fr, ID:12G Length:20 cm

No Star - No Sampling, Check Availability before confirming order of Non Running products.

* - No Sampling Available, Check Availability of Running Products.

** - Sampling Available, No need to check availability of highly moving products, Easy Available.



Double Lumen Curved Catheter

We develop, manufacture, market and support cutting edge vascular access devices and accessories to meet the clinical needs of the medical industry, particularly in the field of interventional cardiology and dialysis. Our company's engineering and applications expertise provides superior products whose progressive designs accommodate advances in medical device and whose quality anticipates the requirements of our professional clients and the patients they serve.

Currently one of the largest manufacturer & importer in the field of Cardiology, Nephrology & Surgical Disposables in INDIA.


Double Lumen Curved Catheter

MECCA HEALTHCARE P LIMITEDs offers Catheter & Kit. It covers almost every size in pediatric and adult length. It comes in single lumen, double lumen & triple lumen with kink resistive catheter.

Standard Kit Contains :

- Guidewire
- Single Lumen Catheter
- Introducer Needle
- Injection Cap
- Dialator

Product Feature



Double D Design

Long term dialysis catheter is made of double D design tubing, which provide excellent flow with low pressure.




Various of Lengths & tip-to-cuff distance

Our long range of catheters help in securing the desired catheter for diverse patient sizes.



Carbothane Material

With second generation Carbothane polyurethane material which facilitates percutaneous placement also this resists collapse during treatment. Which also gives our catheter excellent flexibility.



Infinity tip transition

Our Peel away catheter comes with unique infinity tip transition, which provide excellent insertion and helps in minimizing trauma .



Cuff

Long term dialysis catheter has 9mm of cuff length.



New Tunneler

With atraumatic tip and back end our tunneler helps in minimizing trauma.

Specifications

Double Lumen Curved Catheter


REF	Catheter Length Specification
MI/UDC/DL/01	7 Fr, ID:14G Length:13 cm
MI/UDC/DL/02	7 Fr, ID:14G Length:15 cm
MI/UDC/DL/03	7 Fr, ID:14G Length:16 cm
MI/UDC/DL/04	7 Fr, ID:14G Length:20 cm
MI/UDC/DL/05	8 Fr, ID:12G Length:13 cm
MI/UDC/DL/06	8 Fr, ID:12G Length:15 cm
MI/UDC/DL/07	8 Fr, ID:12G Length:16 cm
MI/UDC/DL/08	8 Fr, ID:12G Length:20 cm

No Star - No Sampling, Check Availability before confirming order of Non Running products.

* No Sampling Available, Check Availability of Running Products.

** Sampling Available, No need to check availability of highly moving products, Easy Available.

Introducer Needle



Introducer Needle

REF	Diameter	Length	Hub Type	Hub Colour
MI/IN/16G70SH	16 Ga	70 mm	Square Hub	Pink
MI/IN/16G70SHW	16 Ga	70 mm	Square Hub	White
MI/IN/16G70R	16 Ga	70 mm	Round Hub	Pink
MI/IN/18G70SH	18 Ga**	70 mm	Square Hub	Pink
MI/IN/18G70SHW	18 Ga**	70 mm	Square Hub	White
MI/IN/18G70R	18 Ga**	70 mm	Round Hub	Pink
MI/IN/18G70Y	18 Ga	70 mm	Y Shape Hub	Blue
MI/IN/20G40SH	20 Ga	40 mm	Square Hub	Pink
MI/IN/20G40SHW	20 Ga*	40 mm	Square Hub	White
MI/IN/20G40Y	20 Ga	40 mm	Y Shape Hub	Blue

No Star - No Sampling, Check Availability before confirming order of Non Running products.

* No Sampling Available, Check Availability of Running Products.

** Sampling Available, No need to check availability of highly moving products, Easy Available.

We develop, manufacture, market and support cutting edge vascular access devices and accessories to meet the clinical needs of the medical industry, particularly in the field of interventional cardiology and dialysis. Our company's engineering and applications expertise provides superior products whose progressive designs accommodate advances in medical device and whose quality anticipates the requirements of our professional clients and the patients they serve.

Introducer Needle

MECCA HEALTHCARE P LIMITED offer a needle with smoother surface, which provide less friction as well as guarantee consistent performance. Our needle has enhanced sharpened tip with polymer coated cannula which provide better and smoother penetration during Seldinger Technique. Flawless metal surface of cannula made of SS304V. We offer Introducer Needle with thin outer wall which help get better inner diameter of cannula to provide more space to introduce MIGwire. MECCA HEALTHCARE P LIMITED offers needle in Square Hub , Round Hub and Y Shape. (This helps to provide the dual port for insertion from single penetration.)

Transducer Protector



Features

- Antibacterial hydrophobic air filter.
- Protective hydrophobic barrier.
- Choice of media to accommodate sterilization methods.
- Pouch Packing.
- Housings designed to withstand pneumatic pulsing.

Transducer Protector

REF	Diameter
MI/TP/01	Standard **

No Star - No Sampling, Check Availability before confirming order of Non Running products.

* No Sampling Available, Check Availability of Running Products.

** Sampling Available, No need to check availability of highly moving products, Easy Available.

We develop manufacture, market and support cutting edge vascular access devices and accessories to meet the clinical needs of the medical industry, particularly in the field of interventional cardiology and dialysis. Our company's engineering and applications expertise provides superior products whose progressive designs accommodate advances in medical device and whose quality anticipates the requirements of our professional clients and the patients they serve.

Currently one of the largest manufacturer & importer in the field of Cardiology, Nephrology & Surgical Disposables in INDIA.

Transducer Protectors

MECCA HEALTHCARE P LIMITED Transducer Protectors are an effective barrier against bacteria and they are used in hemodialysis blood-lines, Our Transducer Protectors are manufactured with fully automatic technology. The membranes are ultrasonically welded and parts are tested with automatic testing equipment in order to guarantee superior quality of the parts. It can be fitted in any blood line or other medical device, covering the market requirement.

A.V. Fistula Needle



Fistula Needle

REF	Diameter	Length	Wing Type	Wing Colour
MI/FNC/RW/16GX1	16 Ga	01 inch	Rotating Wing	Green
MI/FNC/FW/16GX1	16 Ga	01 inch	Fixed Wing	Green
MI/FNC/RW/17GX1	17 Ga	01 inch	Rotating Wing	Orange
MI/FNC/FW/17GX1	17 Ga	01 inch	Fixed Wing	Orange
MI/FNC/RW/15GX1	15 Ga	01 inch	Rotating Wing	Blue
MI/FNC/FW/15GX1	15 Ga	01 inch	Fixed Wing	Blue

No Star - No Sampling, Check Availability before confirming order of Non Running products.

* No Sampling Available, Check Availability of Running Products.

** Sampling Available, No need to check availability of highly moving products, Easy Available.

We develop, manufacture, market and support cutting edge vascular access devices and accessories to meet the clinical needs of the medical industry, particularly in the field of interventional cardiology and dialysis. Our company's engineering and applications expertise provides superior products whose progressive designs accommodate advances in medical device and whose quality anticipates the requirements of our professional clients and the patients they serve.

Currently one of the largest manufacturer & importer in the field of Cardiology, Nephrology & Surgical Disposables in INDIA.

A.V. Fistula Needle

MECCA HEALTHCARE P LIMITED Fistula Needle is used to access the blood vessels for dialysis process. Ultrathin tube and sharp needle reduces the pain on insertion. Soft and elastic tubing is convenient for the patient. Rotatable and fixed wings are easy to handle and control. Screw adapter fit well with blood lines. A.V. Fistula Needle is made to satisfy all clinical requirements.